

3D Files

 **perfect**
www.dental-perfect.com

Extreme Flexibility
Super-elasticity



Make
Endodontics
Easy
with 3D Files

3D Files

**EXTREME
FLEXIBILITY**

**SUPER-
ELASTICITY**

Size	Taper	Torque
● S30	O4	1.0 N·CM
● F25	O2	1.0 N·CM
● F30	O2	1.0 N·CM

Cross section: 

Rotary speed: 800RPM

**REVOLUTIONARY
DESIGN BY
DENTAL PERFECT**

3D FILES adapts perfectly to the shape of the root canal and provide unprecedented cleaning thoroughness

The file processes the remaining uncleaned sections of the channel with scraping movements, and the turbulence it creates increases the efficiency of the irrigant - this allows you to process a large surface area of dentin walls, effectively removing biofilm.

