

Apex Locator  
OLED Screen  
Portable

**perfect**<sup>®</sup>



**ZR-Rap**

ENDO MOTOR



ENDO  
MOTOR

# ZR-Rap



High torque



OLED screen



340°adjustable  
contra angle



1500 mAh  
lithium battery



Low noise

## ◎ Wireless endo motor with built in APEX locator

Wireless Endo motor is generally built with APEX locator and hand-held motor with root-canal length measurement function. It's more convenient to operate relying on multiple functions therein.

## ◎ High efficiency contra angle low power consumption

Relying on higher inverter efficiency and low power consumption, with adoption of high-efficiency curved head and German-imported bearings, it's durable besides of lower power consumption and energy saving.

## ◎ Apex localization activated while preparation

in preparation, it can realize positioning and activation upon top end. In course of various root-canal-therapy preparation stages, including expansion, shape and clean of root canal, it can simultaneously measure length of root canal, which is quite secure.

## ◎ OLED screen that can change its display orientation

With application of digital OLED display screen, it has clear and visual digital displaying effect. Adjustable OLED can also display direction, quite fit for left & right hand synchronized operation.