

Wave Q

 perfect[®]
www.dental-perfect.com

Sonic Irrigation Activator

Wave Q is a specially designed Endo Activator for dentists. With the working tip, the flushing solution is activated through the working tip's vibration to clean the root canal.

13,000
vibrations per minute

62g
holding weight

90 days
battery duration

20mm
holding size



Wave Q

	SIZE	176 X 22 X 22mm
	WEIGHT	60±2G
	BATTERY	3.7V / 800mAh lithium battery
	OUTPUT	5V DC, 2A
	INPUT	AC 100-240~50/60HZ
	MATERIALS	Contra angle: plastic and silicone Working tip: POM plastic

Features

Easy to use

3 different sonic cleaning frequencies.
Easy to install and disassemble in one step.

Avoid breakage

Polymer material(POM) material working tips
avoid breakage during cleaning.

Ergonomic design

Only 62g and the grip part diameter is 20mm
It is the most comfortable size for holding.

Safe and hygienic

Steam sterilizable contra angle avoids cross-infection
Sterilization condition: 135°C, 4 minutes, 0.22Mpa

