

STREAM

Traditional anesthetic procedure is manual controlled and the extrusion speed is unstable, which made the poor absorption effect and the painful feeling of patient.

Stream X is an assistor booster which is specially designed in length and weight, the length is small and the center of gravity of the product is in the front. It relieves patients' pain and doctors' fatigue.

MODE	TIME
Fast Mode	1min
Slow Mode	2min
PDL Mode	4min
Custom Mode	1-4min

Length

166
mm

Weight

72
g

Reflow Time

30
s

Battery Duration

30
days



ASSISTOR
BOOSTER

Eperfect[®]

www.dental-perfect.com

STREAM X

 perfect[®]
www.dental-perfect.com

FEATURES

★ Custom Injection Setting

Stream X is equipped with 4 modes of injection time which varies from 1 to 4 minutes.

★ Short Reflow Setting

The reflow time is just 30 seconds. It's really convenient for the doctor to start the next operation.

★ Comfortable to Hold

Ergonomic & comfortable designed, weight 72g, holding size 20mm, the length is small and the center of gravity is in the front.

★ Durable and Efficient

The 350mAh lithium battery can be continuous used for 30 days(6 root canals per day).

INFORMATION

SIZE	166x21x33mm
WEIGHT	73±2g
BATTERY	3.7V/350mAh lithium battery
BATTERY DURATION	30 days
OUTPUT POWER	5V DC, 2A
INPUT POWER	100-240~50/60Hz
CHARGING TIME	2h(Recommended to charge for more than 3 hours for the first time.)

