

# GP

**perfect**<sup>®</sup>  
www.dental-perfect.com

## Gutta Percha Cutting Pen

The GP Flow is specially designed for dentists to soften or cut the gutta-percha during RCT by heating the working tip.



# FLOW

# GP FLOW

**perfect**<sup>®</sup>  
www.dental-perfect.com

Gutta Percha Cutting Pen

## Features

### Efficient

Rapid heating in 0.2s at 4 levels temperature: 150°C, 180°C, 200°C, 230°C

### Easy to use

The three-button operation logic is simple, making it easy to use.

### Durable

Equipped with 2600mAh large battery and can be continuous used for 30 days

### Comfortable

Ergonomic design: comfortable holding weight 95g and size 20mm, and the working tip is 360° rotatable.

## Information

SIZE	174 x 25 x 25mm
WORKING TIP SIZE	40/.04; 45/.04; 50/.04; 55/.08
HEATING TEMPERATURE	150°C, 180°C, 200°C, 230°C
WEIGHT	95±2g
BATTERY	3.7V / 2600mAh lithium battery
OUTPUT POWER SUPPLY	5V DC, 2A
INPUT POWER SUPPLY	AC 100-240~50/60Hz
STERILIZATION CONDITIONS	135°C, 4minutes, 0.22Mpa

