

# ENDO MASTER

*Brushless Endo Motor*



**E-perfect®**

# FEATURES



## HIGH SPEED BRUSHLESS MOTOR

Using high-performance brushless motors, with high torque, swift response, precise output, and more efficient preparation

Torque Range: 0.4-5.0Ncm  
Speed range: 100-3000rpm



## CUSTOMIZABLE SYSTEM

Customizable preset system file, more convenient to use




## 5.0" IPS TOUCH SCREEN

The base device is equipped with a 5-inch IPS touch screen  
Large screen display, clear at a glance  
Support touch control, intuitive and convenient operation



## LONG ENDURANCE

Handle battery capacity: 1200mAh  
Base battery capacity: 3000mAh



Handle supports wireless charging or Type-c charging



## EASY TO USE



Lighting prompts, human-computer interaction



Audible reminder, Intuitive and clear



Bluetooth connection, wireless operation



# INFORMATION

VELOCITY RANGE	100-3000rpm, brushless motor
TORQUE RANGE	0.4-5.0Ncm
BASE BATTERY CAPACITY	7.4V, 2600mAh
HANDLE BATTERY CAPACITY	3.7V, 1200mAh
ANGLE HANDPIECE	6: 1 Gold and gray
FOLLOW UP SURVEY	Single test and root test in one
HANDLE WEIGHT	About 140g
CHARGING TIME	About 3h

