

Kiddo

Design for children




Scientific operating
sequence
Excellent cutting
efficiency




Kiddo

FEATURES

Material	Niti Gold
Length	16 mm
Sterilization	≤ 134°C
Package	6 pcs per box/blister assorted/ single size

 Special tip design to avoid irritation and damage to the periapical tissue of primary teeth.

 The length is only 16 mm, which is suitable for most children's deciduous teeth root canal length.

INFORMATION

Sample	Size	Length	Torque	Taper	Speed	Color
	#15	16 (mm)	1.50N/cm	0.04	200-350(rpm)	○
	#20	16 (mm)	1.50N/cm	0.04	200-350(rpm)	●
	#25	16 (mm)	1.50N/cm	0.04	200-350(rpm)	●
	#30	16 (mm)	1.50N/cm	0.04	200-350(rpm)	●
	#25	16 (mm)	1.50N/cm	0.04	200-350(rpm)	●
	#30	16 (mm)	1.50N/cm	0.04	200-350(rpm)	●