

 **perfect**®

Wireless 2-in-1 motor

BR
RAP

360°

Rotatable contra-angle

- Industrially Modern
- Exquisitely Made
- Wireless endo Motor



BR-RAP

Wireless 2-in-1 motor

FEATURES



Two-in-one root canal motor

Root canal test alone; It can be shaped while measuring. Two-in-one motor for root testing, precise fit, efficient preparation.



Brushless motors

High-performance brushless motor, low heat generation, low vibration, low noise, smooth operation.



360° rotatable contra-angle motor shank

The contra-angle has a 6:1 transmission ratio and can be freely rotated by 360°. The direction can be adjusted according to the tooth position, which is flexible and versatile.



Wide torque and speed range

Rotational speed: 100-1200rpm.
Easy root canal retreatment, up to 1200 rpm.



• Wireless grip

• 6:1 reduction ratio

• Torque : 0.4-5.0Ncm

• Reciprocating motion, angle adjustable

• Dynamic control, effective prevention of needle breakage

• Large battery capacity, long life without stress