

# Wave Q

 perfect®  
www.dental-perfect.com

## Sonic Irrigation Activator

Wave Q is a specially designed Endo Activator for dentists. With the working tip, the flushing solution is activated through the working tip's vibration to clean the root canal.

**13,000**

vibrations per minute

**62g**

holding weight

**90 days**

battery duration

**20mm**

holding size



# Wave Q

	<b>SIZE</b>	176 X 22 X 22mm
	<b>WEIGHT</b>	60±2G
	<b>BATTERY</b>	3.7V / 800mAh lithium battery
	<b>OUTPUT</b>	5V DC, 2A
	<b>INPUT</b>	AC 100-240~50/60HZ
	<b>MATERIALS</b>	Contra angle: plastic and silicone Working tip: POM plastic

## Features

### Easy to use

3 different sonic cleaning frequencies.  
Easy to install and disassemble in one step.

### Avoid breakage

Polymer material(POM) material working tips  
avoid breakage during cleaning.

### Ergonomic design

Only 62g and the grip part diameter is 20mm  
It is the most comfortable size for holding.

### Safe and hygienic

Steam sterilizable contra angle avoids cross-infection  
Sterilization condition: 135°C, 4 minutes, 0.22Mpa

