

IRRIGATION ACTIVATOR

perfect[®]
www.dental-perfect.com

Ultra Dancer is a specially designed Ultrasonic Endo activate device for dentists. With the working tip, the flushing solution is activated through the working tip's vibration to clean the root canal.

ULTRA DANCER.



✦ 45000 vibrations/s (45KHz) | ✦ 20 days battery duration

©2026 DENTAL PERFECT Co., All Rights Reserved

Ultra Dancer

FEATURES

-  **Efficient**
Ultrasonic frequency: 45KHz
Vibrates 45,000 times in 1s
-  **Easy to use**
3 modes of power · Time customization
Removable and sterilizable handpiece head
-  **Durable**
Equipped with 800mAh large battery and
can be continuous used for 20 days(6 root
canals per day, 1 min per canal)
-  **Comfortable**
Ergonomic & comfortable design
weight 83g, holding size 24mm
-  **Multifunctional**
Suitable for daily irrigation
and additionally supporting
separated file retrieval



INFORMATION

WORKING TIPS SPECIFICATION

| Model | Taper | Working length | Tip diameter |
|-------|-------|----------------|--------------|
| E21 | 0.02 | 21mm | 0.2mm |
| E98 | 0.02 | 18mm | 0.2mm |
| EL96 | 0.015 | 21mm | 0.2mm |

SIZE 192*24*24mm

WEIGHT 83±2g

WORKING TIP MATERIAL Stainless Steel

BATTERY 3.7V/800mAh lithium battery

OUTPUT POWER 5V DC, 2A

INPUT POWER AC100-240~50/60Hz

CHARGING TIME 2.5h(Recommended to charge for more than 3 hours for the first time.)

COMPONENTS

